EN This unique rinsing solution is formulated with OcuPure® preservative to be gentle in your eyes.

LENS PLUS® OcuPure® Saline is the natural choice for rinsing your lenses after disinfection and is also ideal for dissolving protein remover tablets.

Instructions for Use: Always wash and dry your hands thoroughly before handling your contact lenses. After cleaning and disinfecting, rinse your lenses thoroughly with LENS PLUS® OcuPure® Saline prior to insertion. Follow the manufacturer's instructions to use protein remover tablets correctly.

Precautions:

All instructions must be followed. Failure to do so may result in serious eye infection or injury.

- Do not use if tamper-evident seal is broken or missing
- Do not use if allergic to any of the ingredients
- If eye irritation occurs, discontinue lens wear and consult your eye care practitioner
- Always use fresh solution
- Do not use after the expiry date
- Close the bottle tightly immediately after use
- To avoid contamination, do not allow the bottle tip to touch any surface
- Discard any remaining solution within 90 days after first opening the bottle
- Store between 5°C and 25°C
- Keep out of the reach of children
- Do not rinse your lenses or lens case with tap water
- Do not swallow the solution

NOT FOR CONTACT LENS DISINFECTION

Contains:

Sodium chloride, boric acid, OcuPure® preservative (sodium chlorite) - 0.005% and purified water.

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Symbols

Consult instructions for use	Manufacturer	Batch code
Sterile-Aseptic technique	Use by	Temperature limitation
Do not use if package is damaged	Discard any remaining solution 90 days after first opening the bottle	

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■ AMO Ireland
Block B, Liffey Valley Office Campus
Quarryvale, Co Dublin, Ireland

LENS PLUS, LENS PLUS OcuPure and OcuPure are trademarks of Johnson & Johnson Surgical Vision, Inc.

© Johnson & Johnson Surgical Vision, Inc. 2018

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