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Moisturising drops that clean, lubricate and hydrate all types of contact lenses as you blink, with no need to remove them. Moisturising drops with cleaning action.

APPLICATIONS:

Avizor Cleaning Drops are more than just moisturising drops because they hydrate and clean all types of contact lenses while the lenses are worn.

Due to their special formula and cleaning action, they are particularly recommended for users who wear their lenses over long periods of time or who have deposit build-up problems.

INSTRUCTIONS FOR USE:

- Without removing the contact lens, apply one or two drops into each eye and blink several times.
- **Avizor Cleaning Drops** can be used up to four times a day.

Last revision date: 2009/11

PRECAUTIONS:

- Do not use the product after it has been open more than 60 days or after the expiry date.
- Maintain the bottle well closed after use.
- Do not use the product if the safety seal is missing or broken.
- Do not ingest the solution and keep away from children's reach.
- Keep at room temperature.
- Do not allow the dropper to come into contact with any surface to avoid contamination.
- In case of irritation, consult your contact lens specialist.
- Before using any ophthalmic medicine, inform your doctor you are a contact lens user.
- Do not use if you are allergic to any of the components.
- Do not use after the expiration date.

COMPOSITION:

An isotonic, buffered aqueous solution containing Poloxamer, EDTA, Povidone and Polyhexanide 0,0001%.

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