PROS888.2

lacrifresh®
ocu-dry
AVIZOR

0,30%

PRESERVATIVE FREE

Lubricating and hypotonic ophthalmic solution to relieve the symptoms of moderate to severe ocular dryness

HY
Hypotonic solution
Sodium Hyaluronate
0,30% ▲▲▲

gb Lubricating and hypotonic ophthalmic solution to relieve the symptoms of moderate to severe ocular dryness. **PRESERVATIVE FREE. DESCRIPTION: Lacrifresh Ocu-Dry 0.3%** is a lubricating, hypotonic and multidose preservative-free ophthalmic solution with a combination of wetting polymers (sodium hyaluronate, povidone and glycerine) and essential electrolytes, intended to relieve the symptoms of moderate to severe ocular dryness and lubricate all contact lens materials. Due to its gel texture, it is a very suitable product to use also during the night, extending its effect while sleeping. **USAGE: Lacrifresh Ocu-Dry 0.3%** is indicated for: • Prolonged relief of the symptoms of moderate to severe ocular dryness, and other discomforts that are not clinical pathologies, such as a burning, itching, eye strain or the feeling of a foreign body on the surface of the eye, as a consequence of exposure to dusty, smoky, air conditioned or heated environments, prolonged use of computer, hormonal imbalance (menopause), refractive surgery (LASIK), among others. • Lubricating solution for contact lenses. • Due to its composition rich in Sodium Hyaluronate and its high viscosity, it is a very suitable product to be used before going to bed, keeping the eye hydrated throughout the night. **INSTRUCTIONS FOR USE: 1.** Wash, rinse and dry hands before using this product. **2.** Take a bottle of **Lacrifresh Ocu-Dry 0.3%**, pull the safety seal anticlockwise, and remove the protective cap. **3.** Direct use: Apply one or two drops onto each eye and blink several times. **4.** Use with contact lenses: Apply one or two drops on the external surface of each lens during wear or before taken them out, or apply one or two drops into the internal part of each lens before placing on them. **5.** After using the medication, place the cap back on. **6.** Use the product, according to your needs, up to 8 times a day. *Last revision date: 2018/06.* **PRECAUTIONS:** • Do not ingest the solution, and keep out of children's reach. • Do not use the product if you are allergic to any of its ingredients. • Never use tap water to rinse the lens case or the lenses. • If you experience ocular irritation or an adverse reaction, please consult a specialist. • Immediately after instilling it, it is possible to have blurred vision. Wait for a clear vision before performing activities that may involve any risk due to blurred vision. • Once the bottle is opened, use the product within 90 days. • Ensure that the tip of the bottle does not come into contact with the eyes, or any other surface, to avoid contamination. • Do not squeeze the bottle with the protective cap on. • Do not use the product if the safety seal is broken or missing, or the bottle is damaged. • Don't Store the product at room temperature. • Do not use the product after the expiry date. • Keep the bottle tightly closed after each use. • Do not use the product in combination with a pharmacological treatment or any other product that is applied on the ocular surface, without consulting your doctor / eye care practitioner as the functions of both could be altered. **COMPOSITION:** Sodium Hyaluronate 0,30 g, Glycerine, Povidone, Sodium Chloride, Magnesium Chloride, Potassium Chloride, Calcium Chloride, Boric Acid, Borax Decahydrate and Purified Water 100 ml.

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