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CE 0344

3188.20230614.01

ENG Usage instructions
Silicone hydrogel monthly lenses

IMPORTANT: Read this package insert carefully and keep it for later consultation.

INTRODUCTION
This soft lens is made of silicone hydrogel (Definitive 74) with a water content of 74%. Each silicone hydrogel lens is sterile and packaged separately in saline solution. Silicone hydrogel lenses are recommended for daily use for the correction of refractive errors for persons who do not have any diseases of the eye.

Menicon recommends replacing silicone hydrogel lenses after one month. You must follow the cleaning and disinfection plan recommended by the eye care professional or practitioner to clean and disinfect silicone hydrogel lenses after removal. Always use lens care products recommended by your eye care professional or practitioner. Hydrogen peroxide or alcohol-based solutions cannot be used.

CONTRA-INDICATIONS
Do not use silicone hydrogel lenses under the following conditions:
- Acute or sub-acute inflammation or infection of the anterior ocular chamber;
- Any disease, injury or abnormal condition affecting the cornea, conjunctiva or the eyelids;
- Severely insufficient tear production (dry eyes);
- Reduced sensitivity of the cornea;
- Any systemic disease which could negatively affect the eye;
- Allergic reactions that could be triggered by lenses and/or using lens solution;
- Corneal infection;
- Red or irritated eyes;
- Incomplete healing of the cornea.

WARNINGS WHEN WEARING LENSES
Problems resulting from the use of lenses and products for lens care may result in serious eye injury. It is essential that you follow the recommendations of your eye care professional or practitioner and the lens manufacturer (Menicon) regarding the correct use of the lenses and lens care products, including the lens case. Eye problems, including corneal ulcers, may develop very quickly and lead to permanent loss of vision.

You must remove your lenses immediately and contact your eye care professional or practitioner as quickly as possible if your lenses are not comfortable or if you have excessive tear production, changes in your visual acuity or redness in your eyes.

You must remove your lenses immediately and check them carefully if any of the following problems occur:
- Stinging, burning or itchy eyes or other pain in your eyes;
- A persistent scratching feeling;
- Excessive tear production, abnormal ocular secretions, red eyes, reduced vision, blurred vision, flickering, sensitivity to light or dry eyes.

If a lens is damaged, you must not insert it in your eye. Put the lens in the case and immediately contact your eye care professional or practitioner.

If the lens is not damaged, but there is a foreign object on it (eyelash, protein deposit, cosmetics, etc.), you must rinse and disinfect the lens before reinserting it. If the problem persists, remove the lens immediately and contact your eye care professional or practitioner. If any of the problems listed above occur, this may result in a serious situation, such as an infection, a corneal ulcer, neovascularisation or iritis. Remove the lens and immediately contact your eye care professional or practitioner for the correct treatment to prevent serious damage to the eye.

Wearing soft lenses while swimming is not recommended due to the risk of bacterial and/or viral infection.

PREVENTATIVE MEASURES
- Do not touch the lens with your fingernails.
- Thoroughly wash your hands prior to handling the lens. When doing so, avoid using soaps that tend to leave residue behind or towels that have been washed with fabric softener.
- Apply cosmetics only after inserting your lenses. Avoid getting cosmetics, lotions, creams or sprays near your eyes.
- If a lens feels dry in your eye, you should apply a few drops of the by your eye care professional or practitioner recommended lens product, or artificial tears directly into the eye and wait until the lens moves freely on the eye before removing it. If the discomfort lasts longer than 10 minutes, contact your eye care professional or practitioner immediately.
- If the lens dries out during storage, you can no longer wear it.
- Avoid all types of damaging or irritating fumes and smoke while wearing the lenses.
- If chemicals splash into your eye, rinse the eye immediately with tap water. Remove the lenses right away and contact your eye care professional or a doctor.
- Ask your eye care professional or practitioner about wearing lenses during sports or water activities or in certain occupations.
- Inform your doctor that you are a lens wearer.
- Silicone hydrogel lenses should never be worn for longer than one month. After one month, the lens must be disposed of and replaced.
- Never continuously wear the lenses for longer than the period recommended by your eye care professional or practitioner. Silicon hydrogel lenses should never be worn when sleeping.
- Consult your eye care professional or practitioner if your lenses must be stored for longer periods.
- Regular follow-up visits to your eye care professional or practitioner are necessary for all types of lenses to guarantee the health of your eyes.
- Follow all of the instructions in this package insert as well as those of the eye care professional or practitioner.

OPENING THE STERILE PACKAGING
- Silicone hydrogel lenses are steam-sterilised. Never use the lenses if the sterile packaging has been opened or damaged.
- Open the packaging and carefully remove the lens. Inspect the lens to ensure that it is clean and undamaged.
- Rinse the lens with lens care products recommended by your eye care professional or practitioner prior to inserting it for the first time.
- Always start with the same lens to avoid mixing them up.

INSERTING THE LENSES
Menicon recommends the following standard procedure for inserting the lenses:
- Always wash, rinse and dry your hands thoroughly before inserting your lenses.
- Place the lens on the tip of your index finger with the concave side up.
- Make sure the lens is not inside-out by placing it onto your index finger and checking the profile.
- Put a drop of the by your eyecare professional or practitioner recommended lens care product on the inside of the lens.
- Hold your eyelids wide open by pulling up the top lid with the middle finger of your other hand and the lower lid with the middle finger of the hand holding the lens.
- Bring the lens up to your eye while looking straight ahead and place the lens gently onto the eye.
- • Carefully press the lens onto your eye. Release your eyelids and blink gently.
- • Repeat this procedure for the other lens.
- • Empty the lens case, rinse it out with fresh lens care product recommended by your eye care professional or practitioner and close the case until you are ready to use it again.

Note: if you have blurred vision after inserting the lens, check for the following problems:

The lens has shifted:
Look in the mirror to find the lens. With your eye open, carefully move the lens with your finger against the edge of your upper or lower eyelid until the lens is centred.

The lens has not shifted, but one of the following situations has occurred:
- The surface of the lens is dirty;
- The lens is in the wrong eye;
- The lens has been placed onto the eye inside out.

In any of these situations, remove the lens according to the instructions and resolve the issue. If you are unable to resolve the issue, please contact your eye care professional or practitioner.

REMOVING THE LENSES
Menicon recommends the following standard procedure for removing the lenses:
- Always wash, rinse and dry your hands thoroughly before removing your lenses.
- Clean the lens case and fill it with fresh lens care product recommended by your eye care professional or practitioner.
- Always remove the same lens first to avoid mixing them up.
- Look straight ahead in the mirror, hold your head straight and look up until you see the eye in the mirror.
- Pull down the lower lid with the middle finger and place your index finger on the lower edge of the lens.
- Push the lens down to the lower white of your eye.
- Grab the lens carefully using your thumb and index finger.
- Remove the lens.
- Rub the lens clean with lens care product recommended by your eye care professional or practitioner in the palm of your hand for 30 seconds.
- Place the lens in the lens case and ensure that it is fully submerged in the lens care product recommended by your eye care professional or practitioner.
- Repeat the procedure for the other lens.

CARING FOR YOUR LENSES
It is important that you follow the lens care procedure below or the procedure recommended by your eye care professional or practitioner.

When you first get your lenses, you must be able to insert and remove them yourself. You must know and perform the following lens care routine to safely maintain your lenses:
- Always keep the lenses fully submerged in lens care product recommended by your eye care professional or practitioner when you are not wearing them.
- Use lens care product recommended by your eye care professional or practitioner with the lens case provided.
- Replace your lens case on a regular basis. The lens case should be replaced as often as your lens care product.
- Always use fresh lens care product recommended by your eye care professional or practitioner, observe the expiry date(s) provided regardless of whether the product is still sealed or has been opened.
- Do not use tap water, saliva or anything other than the care product recommended by your eye care professional or practitioner for cleaning, rinsing, disinfecting or wetting the lenses/ lens case.

REMINDER: WHEN IN DOUBT, REMOVE THE LENSES AND CONTACT YOUR EYE CARE PROFESSIONAL OR PRACTITIONER.

EXPLANATION OF SYMBOLS:
Batch number / sterilised by steam or dry heat / see instructions for use / expiry date / certification symbol for quality assurance system / manufacturer.

Menicon B.V.
Waanderweg 6
7812 HZ Emmen
The Netherlands

CE 0344

1088.20230614.01

CH REP Company name : Menicon Swiss GmbH (MCH)

Company Address: Dornacherstrasse 27, 4600 Olten, Switzerland

Menicon